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Helping all people
live healthy lives

BD Vacutainer® Order of Draw for Multiple Tube Collections


Designed for Your Safety

Reflects change in NCCLS recommended
Order of Draw (NCCLS H3-A5, Vol 23, No 32, 8.10.2)

* When using a winged blood collection set for venipuncture and a coagulation (citrate) tube is the first specimen tube to be drawn, a discard tube should be drawn first. The discard tube must be used to fill the blood collection set tubing's "dead space" with blood but the discard tube does not need to be completely filled. This important step will ensure maintenance of the proper blood-to-additive ratio of the blood specimen. The discard tube should be a nonadditive or coagulation tube.


Closure Color	Collection Tube	Mix by Inverting
---------------	-----------------	------------------

BD Vacutainer® Blood Collection Tubes (*glass or plastic*)

	• Blood Cultures - SPS	8 to 10 times
	• Citrate Tube*	3 to 4 times
 or 	• BD Vacutainer® SST™ Gel Separator Tube	5 times
	• Serum Tube (<i>glass or plastic</i>)	5 times (plastic) none (glass)
	• Heparin Tube	8 to 10 times
 or 	• BD Vacutainer® PST™ Gel Separator Tube With Heparin	8 to 10 times
 or 	• EDTA Tube	8 to 10 times
	• Fluoride (glucose) Tube	8 to 10 times

**Note: Always follow
your facility's protocol
for order of draw**

Handle all biologic samples and blood collection "sharps" (lancets, needles, luer adapters and blood collection sets) according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since they may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any built-in used needle protector if the blood collection device provides one. BD does not recommend resheathing used needles, but the policies and procedures of your facility may differ and must always be followed. Discard any blood collection "sharps" in biohazard containers approved for their disposal.

BD Global
Technical Services
1.800.631.0174

BD Customer Service
1.888.237.2762
www.bd.com/vacutainer

BD Diagnostics
Preanalytical Systems
1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/vacutainer

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2004 BD.
Printed in USA 06/2004 V55729-4

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1258300672

LabCorp
Laboratory Corporation of America

Shire

Requisition Number:
% r e q n u m b e r %

LABORATORY TEST REQUISITION FORM

Report Test Results To:
Physician:
Site:
LSN:

Contact who completed requisition form:
Name:
Phone # (with country code):

Instructions:
- Complete all information with ball-point pen.
- Submit one legible copy of requisition with each specimen shipment.
- Refer to Project Manual for detailed specimen collection, processing and shipping instructions.

Patient Consent?: ☐ Yes ☐ No
Please note samples may be used for assay qualification.

Physician Name:
Physician Date: DD - MM - YYYY

Tests (Select all that apply)
☐ Anti-Drug Antibody
☐ IgE (for AE only)
☐ TRAP
☐ CHITO
☐ ACE

Visit: Routine - Gaucher

Registry Patient?: ☐ Yes ☐ No If yes, please provide the registry ID.
Registry ID: (Registry patients only) (2 country code / 2 or 3 site # / 4 Patient ID)

Patient ID:

Patient Initials:

Date of Birth: DD - MM - YYYY

Collection Date: DD - MM - YYYY

Collection Time: (24 hour clock)

Adverse Event?: ☐ Yes ☐ No If yes, please describe.
Description of Adv. Event:
Shire Pharmacovigilance or outcome survey groups will be notified of the AE by Labcorp and may follow up with the requesting physician to obtain more information.

Lack of Effect?: ☐ Yes ☐ No

Prior ERT Besides VPRIV?: ☐ Yes ☐ No If yes, please describe.

Date of First VPRIV Infusion: DD - MM - YYYY

For LabCorp Use Only

Expected Container Type	Shipment Condition	Associated Tests	# Containers Expected	# Containers Received	Condition Received		
					Amb	Frz	Refrig
3 mL Cryovial, Serum	Frz	TRAP, CHITO ACE Matrix	1				
3 mL Cryovial, Serum	Frz	ADA, IgE	1				

Accessioner: DE 1: DE 2:

LabCorp Comments:

LabCorp BVBA
Intercity Business Park
Generaal De Wittelaan L11 A2
2800 Mechelen Belgium

LabCorp Central Laboratories (China)
Inc., 203-204 Shidao Soloon
2 Tongji Middle Road, BDA
Beijing 100176 PR China

LabCorp Central Laboratories (Singapore)
Block 33, #03-06, Marsiling Industrial Estate Road 3
Woodlands Central Industrial Estate
Singapore 739256

LabCorp Clinical Trials
750 Walnut Avenue
Cranford, NJ 07016
United States

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2488024525 Laboratory Corporation of America		Requisition Number: % r e q u i s i t i o n n u m b e r %																																		
LABORATORY TEST REQUISITION FORM																																				
Report Test Results To: Physician: _____ Site: _____ LSN: _____ Contact who completed requisition form: Name: _____ Phone # (with country code): _____																																				
Instructions: - Complete all information with ball-point pen. - Submit one legible copy of requisition with each specimen shipment. - Refer to Project Manual for detailed specimen collection, processing and shipping instructions.																																				
Patient Consent?: ☐ Yes ☐ No Please note samples may be used for assay qualification. Physician Name: _____ Physician Date: DD - MM - YY	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;">Visit:</td> <td style="width: 70%;">Routine - Hunter</td> </tr> <tr> <td>Registry Patient?:</td> <td>☐ Yes ☐ No If yes, please provide the registry ID.</td> </tr> <tr> <td>Registry ID: (Registry patients only)</td> <td>DD - MM - YY (2 country code / 3 site # / 4 Patient ID)</td> </tr> <tr> <td>Patient ID:</td> <td>DDMMYY</td> </tr> <tr> <td>Patient Initials:</td> <td>ABC (ABC or A-C)</td> </tr> <tr> <td>Date of Birth:</td> <td>DD - MM - YYYY DD - MM - YYYY</td> </tr> <tr> <td>Collection Date:</td> <td>DD - MM - YYYY DD - MM - YYYY</td> </tr> <tr> <td>Collection Time:</td> <td>HH:MM (24 hour clock)</td> </tr> <tr> <td>Urine Coll. Date:</td> <td>DD - MM - YYYY DD - MM - YYYY</td> </tr> <tr> <td>Urine Coll. Time:</td> <td>HH:MM</td> </tr> <tr> <td>Adverse Event?: ☐ Yes ☐ No If yes, please describe.</td> <td></td> </tr> <tr> <td colspan="2">Description of Adv. Event:</td> </tr> <tr> <td colspan="2">Shire Pharmacovigilance or outcome survey groups will be notified of the AE by Labcorp and may follow up with the requesting physician to obtain more information.</td> </tr> <tr> <td>Lack of Effect?: ☐ Yes ☐ No</td> <td></td> </tr> <tr> <td>Prior ERT: ☐ Yes ☐ No If yes, please describe.</td> <td></td> </tr> <tr> <td>Besides Elaprase?:</td> <td></td> </tr> <tr> <td>Date of First Elaprase Infusion:</td> <td>DD - MM - YYYY DD - MM - YYYY</td> </tr> </table>		Visit:	Routine - Hunter	Registry Patient?:	☐ Yes ☐ No If yes, please provide the registry ID.	Registry ID: (Registry patients only)	DD - MM - YY (2 country code / 3 site # / 4 Patient ID)	Patient ID:	DDMMYY	Patient Initials:	ABC (ABC or A-C)	Date of Birth:	DD - MM - YYYY DD - MM - YYYY	Collection Date:	DD - MM - YYYY DD - MM - YYYY	Collection Time:	HH:MM (24 hour clock)	Urine Coll. Date:	DD - MM - YYYY DD - MM - YYYY	Urine Coll. Time:	HH:MM	Adverse Event?: ☐ Yes ☐ No If yes, please describe.		Description of Adv. Event:		Shire Pharmacovigilance or outcome survey groups will be notified of the AE by Labcorp and may follow up with the requesting physician to obtain more information.		Lack of Effect?: ☐ Yes ☐ No		Prior ERT: ☐ Yes ☐ No If yes, please describe.		Besides Elaprase?:		Date of First Elaprase Infusion:	DD - MM - YYYY DD - MM - YYYY
Visit:	Routine - Hunter																																			
Registry Patient?:	☐ Yes ☐ No If yes, please provide the registry ID.																																			
Registry ID: (Registry patients only)	DD - MM - YY (2 country code / 3 site # / 4 Patient ID)																																			
Patient ID:	DDMMYY																																			
Patient Initials:	ABC (ABC or A-C)																																			
Date of Birth:	DD - MM - YYYY DD - MM - YYYY																																			
Collection Date:	DD - MM - YYYY DD - MM - YYYY																																			
Collection Time:	HH:MM (24 hour clock)																																			
Urine Coll. Date:	DD - MM - YYYY DD - MM - YYYY																																			
Urine Coll. Time:	HH:MM																																			
Adverse Event?: ☐ Yes ☐ No If yes, please describe.																																				
Description of Adv. Event:																																				
Shire Pharmacovigilance or outcome survey groups will be notified of the AE by Labcorp and may follow up with the requesting physician to obtain more information.																																				
Lack of Effect?: ☐ Yes ☐ No																																				
Prior ERT: ☐ Yes ☐ No If yes, please describe.																																				
Besides Elaprase?:																																				
Date of First Elaprase Infusion:	DD - MM - YYYY DD - MM - YYYY																																			
Tests (Select all that apply) ☐ Anti-Drug Antibody ☐ IgE (for AE only) ☐ GAG ☐ Creatinine (for GAG)																																				
For LabCorp Use Only																																				
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Expected Container Type</th> <th>Shipment Condition</th> <th>Associated Tests</th> <th># Containers Expected</th> <th># Containers Received</th> <th>Condition Received</th> </tr> </thead> <tbody> <tr> <td>3 mL Cryovial, Serum</td> <td>Frz</td> <td>ADA, IgE</td> <td>1</td> <td></td> <td></td> </tr> <tr> <td>3 mL Cryovial, Urine</td> <td>Frz</td> <td>GAG</td> <td>1</td> <td></td> <td></td> </tr> <tr> <td>3 mL Cryovial, Urine</td> <td>Frz</td> <td>Creatinine</td> <td>1</td> <td></td> <td></td> </tr> </tbody> </table>	Expected Container Type	Shipment Condition	Associated Tests	# Containers Expected	# Containers Received	Condition Received	3 mL Cryovial, Serum	Frz	ADA, IgE	1			3 mL Cryovial, Urine	Frz	GAG	1			3 mL Cryovial, Urine	Frz	Creatinine	1			<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Accessioner:</th> <th>DE 1:</th> <th>DE 2:</th> </tr> </thead> <tbody> <tr> <td colspan="3" style="height: 40px; vertical-align: bottom;"> LabCorp Comments: </td> </tr> </tbody> </table>		Accessioner:	DE 1:	DE 2:	LabCorp Comments:						
Expected Container Type	Shipment Condition	Associated Tests	# Containers Expected	# Containers Received	Condition Received																															
3 mL Cryovial, Serum	Frz	ADA, IgE	1																																	
3 mL Cryovial, Urine	Frz	GAG	1																																	
3 mL Cryovial, Urine	Frz	Creatinine	1																																	
Accessioner:	DE 1:	DE 2:																																		
LabCorp Comments:																																				
<table style="width: 100%;"> <tr> <td style="width: 33%;"> LabCorp BVBA Intercity Business Park Generaal De Wittelaan L11 A2 2800 Mechelen Belgium </td> <td style="width: 33%;"> LabCorp Central Laboratories (China) Inc., 203-204 Shidao Solaon 2 Tongji Middle Road, BDA Beijing 100176 PR China </td> <td style="width: 33%;"> LabCorp Central Laboratories (Singapore) Pte. Ltd. Block 33, #03-06, Marsiling Industrial Estate Road 3 Woodlands Central Industrial Estate Singapore 739256 </td> </tr> <tr> <td> LabCorp Clinical Trials 750 Walnut Avenue Cranford, NJ 07016 United States </td> <td colspan="2"></td> </tr> </table>			LabCorp BVBA Intercity Business Park Generaal De Wittelaan L11 A2 2800 Mechelen Belgium	LabCorp Central Laboratories (China) Inc., 203-204 Shidao Solaon 2 Tongji Middle Road, BDA Beijing 100176 PR China	LabCorp Central Laboratories (Singapore) Pte. Ltd. Block 33, #03-06, Marsiling Industrial Estate Road 3 Woodlands Central Industrial Estate Singapore 739256	LabCorp Clinical Trials 750 Walnut Avenue Cranford, NJ 07016 United States																														
LabCorp BVBA Intercity Business Park Generaal De Wittelaan L11 A2 2800 Mechelen Belgium	LabCorp Central Laboratories (China) Inc., 203-204 Shidao Solaon 2 Tongji Middle Road, BDA Beijing 100176 PR China	LabCorp Central Laboratories (Singapore) Pte. Ltd. Block 33, #03-06, Marsiling Industrial Estate Road 3 Woodlands Central Industrial Estate Singapore 739256																																		
LabCorp Clinical Trials 750 Walnut Avenue Cranford, NJ 07016 United States																																				

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Site: _____

LSN: _____

Contact who completed requisition form:

Name: _____

Phone # (with country code): _____

Visit:	Routine - Fabry
Registry Patient?:	☐ Yes ☐ No If yes, please provide the registry ID.
Registry ID:	[] [] - [] [] - [] [] [] [] (Registry patients only) (2 country code / 2 or 3 site # / 4 Patient ID)
Patient ID:	[] [] [] [] [] [] [] []
Patient Initials:	[] [] (ABC or A-C)
Date of Birth:	[] [] - [] [] [] - [] [] [] [] D D - M M M - Y Y Y Y
Collection Date:	[] [] - [] [] [] - [] [] [] [] D D - M M M - Y Y Y Y
Collection Time:	[] [] : [] [] (24 hour clock)
Urine Coll. Date:	[] [] - [] [] [] - [] [] [] [] D D - M M M - Y Y Y Y
Urine Coll. Time:	[] [] : [] []
Adverse Event?:	☐ Yes ☐ No If yes, please describe.
Description of Adv. Event:	
Shire Pharmacovigilance or outcome survey groups will be notified of the AE by Labcorp and may follow up with the requesting physician to obtain more information.	
Lack of Effect?:	☐ Yes ☐ No
Prior ERT Besides Replagal:	☐ Yes ☐ No If yes, please describe.
Date of First Replagal Infusion:	[] [] - [] [] [] - [] [] [] [] D D - M M M - Y Y Y Y

Instructions:

- Complete all information with ball-point pen.
- Submit one legible copy of requisition with each specimen shipment.
- Refer to Project Manual for detailed specimen collection, processing and shipping instructions.

Patient Consent?: ☐ Yes ☐ No

Please note samples may be used for assay qualification.

Physician Name: _____

Physician Date: _____

Tests (Select all that apply)

☐ Anti-Drug Antibody

☐ IgE (for AE only)

☐ Creatinine (for urine GB-3)

☐ GB-3 (urine)

☐ GB-3/Lyso GB-3 (plasma)

For LabCorp Use Only

Expected Container Type	Shipment Condition	Associated Tests	# Containers Expected	# Containers Received	Condition Received		
					Amb	Frz	Refrig
3 mL Cryovial, Serum	Frz	Anti-Drug Ab, IgE	1				
10 mL Cryovial, Urine	Frz	GB-3	1				
3 mL Cryovial, Urine	Frz	Creatinine	1				
3 mL Cryovial, Plasma	Frz	GB-3/Lyso GB-3	1				

Accessioner: _____

LabCorp Comments: _____

DE 1: _____

DE 2: _____

LabCorp BVBA
 Intercity Business Park
 Generaal De Wittelaan L11 A2
 2800 Mechelen Belgium

LabCorp Central Laboratories (China)
 Inc., 203-204 Shidao Soloon
 2 Tongji Middle Road, BDA
 Beijing 100176 PR China

LabCorp Central Laboratories (Singapore)
 Block 33, #03-06, Marsiling Industrial Estate Road 3
 Woodlands Central Industrial Estate
 Singapore 739256

LabCorp Clinical Trials
 750 Walnut Avenue
 Cranford, NJ 07016
 United States

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REQUEST FOR DRY ICE SHIPPING BOX DELIVERY (FAX LATEST BY 14:00 C.E.T. 1 DAY PRIOR TO REQUESTED DELIVERY DATE)			
Please fax this form to +32 (0)2 307 9150 or e-mail to brulosond@dhl.com You can contact DHL at +32 2 307 9138 For contact details of the Esoterix monitor team consult your manual			
Reference for billing (Esoterix GCP/ESN)			
Protocol number			
Site number (= Esoterix Account Number)			
Investigator Name			
MARK YOUR REQUIREMENT ON AMOUNT OF SAMPLES TO RETURN: 1 - 32 SAMPLES ☐ (1 small box (TYPE 8.1)) 33 - 64 SAMPLES ☐ (1 big box (TYPE 13.1))		PRE12 REQUIRED ? ☐ Sponsor approval required	
If you ship more than 64 samples, please make a combination of both box types based on the total number of samples to ship			
REQUESTED DELIVERY DATE: YOUR SITE HAS A 24 hrs SERVICE ONLY USE FOR A SINGLE DELIVERY DATE - IN MULTIPLE DELIVERIES USE SEPARATE FORMS PER DELIVERY DATE			
DRY ICE DELIVERY POSSIBLE ON: TUESDAY WEDNESDAY THURSDAY	REQUESTED DELIVERY DATE DD / MM / YY	DHL INTERNAL USE: BOX TYPE 8.1 / 13.1 (Sent from dispatch area on Monday) (Sent from dispatch area on Tuesday) (Sent from dispatch area on Wednesday)	
Hospital name			
Department name			
Attention of			
Address			
City		Post/Zip Code:	
Country			
Phone			
Fax			
Email			
Date	Ordered by		
IF YOU SHIP SAMPLES ON FRIDAY DO NOT FORGET YOUR SATURDAY DELIVERY STICKER			
To be completed by DHL (for internal use): Dispatch AWB Return AWB Delivery Station			

Dry Ice Box 24 hrs zone - solution 8.1 / solution 13.1

CT-BG-KL-FCD-19 Revision 1

TNT Express - Dry Ice order form MEDPAK - 24 hrs site			
Please FAX this form to +32 (0) 2 754.49.47 or e-mail to Be.Dryice@tnt.com at the latest by 14:00 C.E.T. 1 day BEFORE the desired delivery date			
If you have any problems receiving the dry ice please contact your Esoterix Monitor team IMMEDIATELY (contact details see manual).			
Please complete all the information below for delivery of Medpak & Dry Ice by TNT			
Name of Central Laboratory		Esoterix Clinical Trials Services	
Protocol Number + (GCP / ESN Code)		TNT account BE 67675	
Site Number			
Name Main Investigator			
Hospital/Clinic			
Department/Wing			
Street Name and nr (Do not use PO Box Numbers)			
Postal code + City			
Country			
Name Contact Person for the shipment			
Contact Telephone Number			
Fax Number (Fill out for TNT confirmation)			
Ordered By			
Comments			
Dry Ice and packaging required at site on: Date: _____ Day: _____		Clearly fill in the number of boxes needed: ☐ 1 - 12 samples = 1 box MedPak 5 ☐ 12 - 24 samples = 1 box MedPak 10	
PRE12 REQUIRED ? ☐ Sponsor approval required		IF YOU SHIP SAMPLES ON FRIDAY DO NOT FORGET YOUR SATURDAY DELIVERY STICKER	
Following boxes for TNT office use only			
Order received by			
Time received (fax time)			
Dry Ice form 24 hrs zone			
CT-BG-KL-FCD-20 Revision 1			

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The image shows a DHL Shipment Air Waybill form. Key sections include:

- Payer account number and insurance details:** Fields for Payer Account No., Shipper's account number, and insurance options.
- Shipper details:** Fields for Shipper's name, address, and contact information.
- Destination code:** A field for the destination code.
- Shipment details:** Fields for Total number of packages, Total weight, and dimensions.
- Full description of contents:** A field for describing the contents of the shipment.
- Signature:** A field for the shipper's signature, circled in blue with an arrow pointing to the instruction 'Datez et signez'.

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The image shows a TNT Shipment Air Waybill form. Key sections include:

- From (Collection Address):** Fields for Name, Collection Address, City, Province, Postal/Zip Code, Country, Contact Name, and Tel No.
- To (Receiver Address):** Fields for Name, Receiver Address, City, Province, Postal/Zip Code, Country, Contact Name, and Tel No.
- Signature:** A field for the shipper's signature, circled in blue with an arrow pointing to the instruction 'Datez et signez'.

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Non applicable

LABORATORY ACCREDITATION

A valid Laboratory Accreditation is provided with the initial supply.

For online accreditation information, access:

http://www.labcorpclinicaltrials.com/global_locations_and_accreditations

LABORATORY ACCREDITATION – MECHELEN

For more information about the accreditation and a current list of accredited tests please refer to the following website:
www.belac.be.

 Organisme belge d'Accréditation Belgische Accreditatieinstelling Belgische Akkreditierungsstelle Belgian Accreditation Body Signatory to EA, ILAC and JAF Multilateral Agreements	Accreditation Certificate No. 299-TEST In compliance with the provisions of the Royal Decree of 31 January 2006 setting up BELAC, the Accreditation Board hereby declares, that the test laboratory LABCORP CLINICAL TRIALS INTERCITY BUSINESS PARK Generaal De Wittelaan, L11 A2 2800 MECHELEN - Belgium has the competence to perform the tests as described in the annex which is an integral part of the present certificate, in accordance with the requirements of the standard NBN EN ISO/IEC 17025:2005. The present accreditation is the subject of regular surveillance in order to confirm the compliance with the accreditation conditions. The Chair of the Accreditation Board BELAC,  Nicole MEURÉE-VANLAETHEM Issue date : 2014-06-17 Validity date : 2014-12-17 Original version of this certificate is in Dutch. 
---	---

Advancing Excellence

Accredited
Laboratory

The College of American Pathologists

certifies that the laboratory named below

**LabCorp Clinical Trials
Laboratory
Mechelen, Belgium
Raf Schepers, PhD**

LAP Number: 6967001
AU-ID: 1197429

*has met all applicable standards for accreditation and
is hereby accredited by the College of American Pathologists'
Laboratory Accreditation Program. Reinspection should occur prior
to January 30, 2015 to maintain accreditation.*

Accreditation does not automatically survive a change in director, ownership,
or location and assumes that all interim requirements are met.

Chair, Commission on Laboratory Accreditation

President, College of American Pathologists


LABORATORY ACCREDITATION

A valid Laboratory Accreditation is provided with the initial supply.

For online accreditation information, access:

http://www.labcorpclinicaltrials.com/global_locations_and_accreditations

LABORATORY ACCREDITATION – CRANFORD, NJ

**The College of American Pathologists**

certifies that the laboratory named below

LabCorp Clinical Trials Services

Laboratory

Cranford, New Jersey

Steven M. Diamond, DO

LAP Number: 7179720

AU-ID: 1387682

CLIA Number: 31D1004205

has met all applicable standards for accreditation and is hereby accredited by the College of American Pathologists' Laboratory Accreditation Program. Reinspection should occur prior to April 15, 2015 to maintain accreditation.

Accreditation does not automatically survive a change in director, ownership, or location and assumes that all interim requirements are met.

Chair, Commission on Laboratory Accreditation

President, College of American Pathologists

CENTERS FOR MEDICARE & MEDICAID SERVICES
CLINICAL LABORATORY IMPROVEMENT AMENDMENTS
CERTIFICATE OF ACCREDITATION


LABORATORY NAME AND ADDRESS LABCORP CLINICAL TRIALS 750 WALNUT AVENUE CRANFORD, NJ 07016 LABORATORY DIRECTOR STEVEN M DIAMOND DO	CLIA ID NUMBER 31D1004205 EFFECTIVE DATE 07/01/2013 EXPIRATION DATE 06/30/2015
--	---

Pursuant to Section 353 of the Public Health Services Act (42 U.S.C. 263a) as revised by the Clinical Laboratory Improvement Amendments (CLIA), the above named laboratory located at the address shown herein (and other approved locations) may accept human specimens for the purposes of performing laboratory examinations or procedures.

This certificate shall be valid until the expiration date above, but is subject to revocation, suspension, limitation, or other sanctions for violation of the Act or the regulations promulgated thereunder.

CENTERS FOR MEDICARE & MEDICAID SERVICES

 Judith A. Yost, Director
 Division of Laboratory Services
 Survey and Certification Group
 Center for Medicaid and State Operations

21 Cmr2_060513

If you currently hold a Certificate of Compliance or Certificate of Accreditation, below is a list of the laboratory specialties/subspecialties you are certified to perform and their effective dates:

LAB CERTIFICATION (CODE)	EFFECTIVE DATE	LAB CERTIFICATION (CODE)	EFFECTIVE DATE
VIROLOGY (140)	04/03/2007		
GENERAL IMMUNOLOGY (220)	07/01/2003		
ROUTINE CHEMISTRY (310)	07/01/2003		
URINALYSIS (320)	07/01/2003		
ENDOCRINOLOGY (330)	07/01/2003		
TOXICOLOGY (340)	08/23/2004		
HEMATOLOGY (400)	07/01/2003		

FOR MORE INFORMATION ABOUT CLIA, VISIT OUR WEBSITE AT WWW.CMS.HHS.GOV/CLIA
 OR CONTACT YOUR LOCAL STATE AGENCY. PLEASE SEE THE REVERSE FOR
 YOUR STATE AGENCY'S ADDRESS AND PHONE NUMBER.
 PLEASE CONTACT YOUR STATE AGENCY FOR ANY CHANGES TO YOUR CURRENT CERTIFICATE.

LABORATORY ACCREDITATION – TANDEM LABS RTP

CENTERS FOR MEDICARE & MEDICAID SERVICES CLINICAL LABORATORY IMPROVEMENT AMENDMENTS CERTIFICATE OF REGISTRATION	
LABORATORY NAME AND ADDRESS	CLIA ID NUMBER
TANDEM LABS - RTP A LABORATORY CORPORATION OF AMERICA HO 2202 ELLIS ROAD, SUITES B, C, D DURHAM, NC 27703-5518 LABORATORY DIRECTOR	34D2074059 EFFECTIVE DATE 02/28/2014 EXPIRATION DATE 02/28/2016
PERCY E JONES M.D.	
Pursuant to Section 355 of the Public Health Services Act (42 U.S.C. 2653a) as revised by the Clinical Laboratory Improvement Amendments (CLIA), the above named laboratory located at the address shown herein (and other approved locations) may accept human specimens for the purposes of performing laboratory examinations or procedures. This certificate shall be valid until the expiration date above, but is subject to revocation, suspension, limitation, or other sanctions for violation of the Act or the regulations promulgated thereunder.	
	 Judith A. Yost, Director Division of Laboratory Services Survey and Certification Group Center for Clinical Standards and Quality

LABORATORY ACCREDITATION – TANDEM LABS NEW JERSEY

CENTERS FOR MEDICARE & MEDICAID SERVICES CLINICAL LABORATORY IMPROVEMENT AMENDMENTS CERTIFICATE OF COMPLIANCE	
LABORATORY NAME AND ADDRESS	CLIA ID NUMBER
TANDEM LABS A DIVISION OF NWT INC 115 SILVIA STREET WEST TRENTON, NJ 08628	31D1067412
LABORATORY DIRECTOR	EFFECTIVE DATE
STEVEN M DIAMOND DO	06/06/2014
	EXPIRATION DATE
	06/05/2016

Pursuant to Section 353 of the Public Health Services Act (42 U.S.C. 263a) as revised by the Clinical Laboratory Improvement Amendments (CLIA), the above named laboratory located at the address shown herein (and other approved locations) may accept human specimens for the purposes of performing laboratory examinations or procedures.

This certificate shall be valid until the expiration date above, but is subject to revocation, suspension, limitation, or other sanctions for violation of the Act or the regulations promulgated thereunder.

 *Judith A. Yost*
Judith A. Yost, Director
Division of Laboratory Services
Survey and Certification Group
Center for Clinical Standards and Quality

44 Cert02_000014

If you currently hold a Certificate of Compliance or Certificate of Accreditation, below is a list of the laboratory specialties/subspecialties you are certified to perform and their effective date:

LAB CERTIFICATION (CODE)	EFFECTIVE DATE	LAB CERTIFICATION (CODE)	EFFECTIVE DATE
GENERAL IMMUNOLOGY (220)	06/06/2008		
TOXICOLOGY (340)	06/06/2008		

FOR MORE INFORMATION ABOUT CLIA, VISIT OUR WEBSITE AT WWW.CMS.GOV/CLIA
OR CONTACT YOUR LOCAL STATE AGENCY. PLEASE SEE THE REVERSE FOR
YOUR STATE AGENCY'S ADDRESS AND PHONE NUMBER.
PLEASE CONTACT YOUR STATE AGENCY FOR ANY CHANGES TO YOUR CURRENT CERTIFICATE.

New Jersey Department of Health
DIVISION OF PUBLIC HEALTH AND ENVIRONMENTAL LABORATORIES

CLINICAL LABORATORY LICENSE

No. 00031016

The following, pursuant to Chapter 166, P.L. of 1975, is hereby authorized to perform the below indicated services:

TANDEM LABS -
115 SILVIA ST
WEST TRENTON, NJ 08638

CUS ID: 0002750
Effective: 01/01/2014
To: 12/31/2014

AUTHORIZED SERVICES		
☐ Urinalysis	☐ Mycology	☐ Chemistry
☐ Bacteriology	☐ Class I ☐ Class III	☐ Limited
☐ Limited	☐ Class II ☐ Class IV	
	☐ Virology	☐ Endocrinology
☐ Mycobacteriology	☒ Diagnostic Immunology	☒ Toxicology
☐ Class I	☐ Syphilis Serology	☐ Cytology
☐ Class II	☒ General Immunology	☐ Collection Station Only
☐ Class III	☐ Hematology	☐ Cytogenetics and/or Tissue Typing
☐ Class IV	☐ Limited	☐ Collection Station Performing Waived Tests Only
☐ Parasitology		☐ Other
☐ Limited	☐ Immunohematology	☐ Limited
	☐ Group and Type Only	

TO BE CONSPICUOUSLY DISPLAYED AT LABORATORY

COMMISSIONER OF HEALTH

CV MEDICAL DIRECTOR - MECHELEN

Raf Schepers

Kapellelaan 89
1860 Meise
Belgium
Nationality: Belgian

April 1st, 2014

Cell: + 32 475 32 20 70
rscheper@gmail.com

EDUCATION

Pharmacy degree, 1995, Katholieke Universiteit Leuven (KUL), Leuven, Belgium
Master of Industrial Pharmaceutics (cum laude), 1996 - KUL, Leuven, Belgium
Clinical Pathology degree (cum maxima laude), 2001, Université Catholique de Louvain (UCL), Brussels, Belgium
Radiation Protection and Radioimmunoassay degree, 2001, UCL, Brussel – Autorisation to use radioisotopes Federal Agency for Nuclear Control, 2010.
Ph.D. Neuroscience, 2007, National Institutes of Health, Baltimore, Maryland, USA, University of Maryland, Baltimore, Maryland, USA
QA degree, 2010, KUL, Leuven, Belgium

EXPERIENCE

LabCorp Clinical Trials – Jan 2014 – current
Global Director LabCorp Clinical Trial laboratories
Senior Laboratory Director LabCorp Clinical Trials Mechelen, Belgium

Centraal Laboratorium Hasselt – Centrum Medische analyses, May 2011 – Dec 2013, Hasselt, Belgium
Clinical Pathologist
Head laboratory operations Toxicology (LCMSMS, GCMS, ICPMS ...)

Medisch Laboratorium Medina, May 2009 – March 2014, Dendermonde and Aalter, Belgium
QA Manager
Clinical Pathologist

BARC N.V., Sept 2008 – April 2009, Bioanalytical Research Company, Zwijnaarde, Belgium
Head laboratory operations RIA and Toxicology
Clinical Pathologist

Post-doctoral Fellow, Aug 2007- Aug 2008, Department of Neurosurgery, **Johns Hopkins University**, Baltimore, Maryland, USA., Principal investigators: Dr. M. Ringkamp, Dr. RA Meyer, Dr. J. Campbell
Post-doctoral research with regards to the discovery of mechanisms of itch and neuropathic pain (in vivo electrophysiology studies in monkey).

Doctoral Fellow, Nov 2003 – Juli 2007, Integrative Neurosciences Branch, **NIDA**, NIH, Baltimore, Maryland, USA, Advisor: Dr. Toni Schippenberg Thesis: Mechanisms of brainstem mu and kappa opioid mediated descending modulation of inflammatory and neuropathic pain. (analytical skills mastered: neurosurgery, brain microdialysis, microbore HPLC with laser fluorescence detection of neurotransmitters)

NIH Exchange Visitor, Oct 2001 – Oct 2003, Chemistry and Drug Metabolism Section, Clinical Pharmacology and Therapeutics Research Branch, **NIDA**, NIH, Baltimore, Maryland, USA, Advisor: Dr. Marilyn Huestis
Clinical researcher: principal investigator of PD and PK studies pertaining various drugs of abuse. (analytical skills mastered: GC-MS / LCMSMS / MALDI-TOF)

RAF SCHEPERS
p1/3

HONOURS

Reviewer: The International Journal of Neuropsychopharmacology
Clinical Pharmacology and Therapeutics
Pharmacology Biochemistry and Behavior

IASP travel award IASP meeting 2008
Democratically elected student representative, Faculty Farmaceutical sciences, KUL, Leuven, Belgium, 1993-1995.

LANGUAGE SKILLS

- trilingual Dutch, English, French
- German: good understanding, basic speaking and writing

PUBLICATIONS

Zapata A, **Schepers RJ**, Pontis S, Gehrke B, Minney V, Cai H, Figueiró-Silva J, Podlesniy P, Trullas R, and Shippenberg TS. Neuronal Pentraxin in the Rostral Ventromedial Medulla contributes to the descending inhibition of tissue injury evoked pain in rodents. Submitted to Journal of Pain

Zapata A, Pontis S, **Schepers RJ**, Wang R, Oh E, Stein A, Bäckman CM, Worley P, Enguita M, Abad MA, Trullas R, Shippenberg TS. Alleviation of neuropathic pain hypersensitivity by inhibiting neuronal pentraxin 1 in the rostral ventromedial medulla. J Neurosci. 2012 Sep 5;32(36):12431-6.

Ringkamp M, **Schepers RJ**, Shimada SG, Johanek LM, Hartke TV, Borzan J, Shim B, LaMotte RH, Meyer RA. A role for nociceptive, myelinated nerve fibers in itch sensation. J Neurosci. 2011 Oct 19;31(42):14841-9.

Ringkamp M, Johanek LM, Borzan J, Hartke TV, Wu G, Pogatzki-Zahn EM, Campbell JN, Shim B, **Schepers RJ**, Meyer RA. PLoS One. Conduction properties distinguish unmyelinated sympathetic efferent fibers and unmyelinated primary afferent fibers in the monkey. 2010 Feb 5;5(2):e9076.

Schepers RJ and Ringkamp M. Thermoreceptors and thermosensitive afferents. Neurosci Biobehav Rev. 2008 Aug 7.

Schepers RJ, Mahoney JL, Zapata A, Chefer V, Shippenberg TS. The effects of local perfusion of DAMGO on extracellular GABA and glutamate concentrations in the rostral ventromedial medulla. J Neurochem. 2008 Feb;104(3):806-17.

Schepers RJ, Mahoney JL, Gehrke BJ, Shippenberg TS. Endogenous Kappa-Opioid Receptor Systems Inhibit Hyperalgesia Associated With Localized Peripheral Inflammation. Pain. 2008; Pain. 2008 Aug 31;138(2):423-39.

Schepers RJ, Mahoney JL, Shippenberg TS. Inflammation-induced changes in rostral ventromedial medulla mu and kappa opioid receptor mediated antinociception. Pain. 2008 Jun;136(3):320-30.

Schepers RJ, Oyler JM, Joseph RE Jr, Cone EJ, Moolchan ET, Huestis MA. Methamphetamine and amphetamine pharmacokinetics in oral fluid and plasma after controlled oral methamphetamine administration to human volunteers. Clin Chem. 2003;49(1):121-32.

Kim I, Barnes AJ, Oyler JM, **Schepers R**, Joseph RE Jr, Cone EJ, Lafko D, Moolchan ET, Huestis MA. Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration. Clin Chem. 2002;48(9):1486-96.

RAE
SCHEPERS
p213

Schepers R, Mourad M, Wagner C, Chaib-Eddour D, Malaise J, Squifflet JP and Wallemacq PE. Evaluation of cystatin C and serum amyloid A (SAA) in the follow up of 22 kidney transplant patients. *Laboratorine Medicina* 2001;10:10-16.

Mourad M, Malaise J, Chaib Eddour D, De Meyer M, Konig J, **Schepers R**, Squifflet JP, Wallemacq P. Pharmacokinetic basis for the efficient and safe use of low-dose mycophenolate mofetil in combination with tacrolimus in kidney transplantation. *Clin Chem.* 2001;47(7):1241-8.

Mourad M, Malaise J, Chaib Eddour D, De Meyer M, Konig J, **Schepers R**, Squifflet JP, Wallemacq P. Correlation of mycophenolic acid pharmacokinetic parameters with side effects in kidney transplant patients treated with mycophenolate mofetil. *Clin Chem.* 2001;47(1):88-94.

Adams E, **Schepers R**, Gathu LW, Kibaya R, Roets E, Hoogmartens J. Liquid chromatographic analysis of a formulation containing polymyxin, gramicidin and neomycin. *Journal of Pharmaceutical and Biomedical Analysis* 1997;15(4):505-511.

Adams E, **Schepers R**, Roets E, Hoogmartens J. Determination of neomycin sulfate by liquid chromatography with pulsed electrochemical detection. *Journal of Chromatography A* 1996;741(2):233-240.

POSTERS AND PRESENTATIONS

Oral presentation, Matrix Assisted Laser Desorption Ionization – Time of Flight Mass Spectrometry: methods and applications for forensic and toxicologic research, Society of Forensic Toxicology (SOFT) meeting, Milwaukee, IL, 2000.

Oral presentation, Methamphetamine and amphetamine pharmacokinetics in oral fluid and plasma after controlled oral methamphetamine administration to human volunteers. SOFT, New Orleans, LA, 2001.

Poster, Inflammation-induced changes in rostral ventromedial medulla mu and kappa opioid receptor mediated antinociception, International Association for the study of Pain (IASP), Sydney, Australia, 2004

Poster, The effects of local perfusion of DAMGO on extracellular GABA and glutamate concentrations in the rostral ventromedial medulla, Society for Neuroscience (SfN), San Diego, USA 2007

Oral presentation, Kappa opioid receptor antinociception: recent advances in basic and clinical research. Blaustein seminars, JHU, 2008

Poster, A subpopulation of mechanosensitive A-fiber nociceptors in monkey respond to pruritic stimuli, IASP, 2008

EXTRACURRICULAR ACTIVITIES

- Organ, first price interdisciplinary competition of music instruments, awarded by the city of Hasselt, 1990
- Bassoon, level 4, Eisden music school, 1990
- Piano, level 4, Hasselt music conservatory, 1993
- Accompanist, choir "Virga Jesse", 1988-1996, quintet Cantabile 1993-1995

RAF SCHEPERS

p 313

CV MEDICAL DIRECTOR – CRANFORD

CURRICULUM VITAE

Steven Diamond, DO, MBA, FCAP

<u>Board Certification:</u>	American Board of Pathology Anatomic and Clinical Pathology (Recertified) Hematology	11/75 1/09 5/77
<u>Board Registry:</u>	America Board of Medicolegal Death Investigators	5/20/02
<u>Education:</u>	Chicago College of Osteopathic Medicine Midwestern University Chicago, IL	DO 6/70
	University of Wisconsin Madison, WI	Path Resident 1970 – 1971
	United States Navy Philadelphia, PA	Pathologist 1971 – 1973
	Thomas Jefferson University Philadelphia, PA	Path Resident 1973 – 1975
	Wharton School of Business University of Pennsylvania Philadelphia, PA	MBA 1992 – 1994
<u>State Licenses:</u>	New Jersey Pennsylvania Missouri New York North Carolina South Carolina Virginia Georgia Ohio Indiana	
<u>Medical Societies:</u>	College of American Pathologists Pennsylvania Medical Society National Association of Medical Examiners American Board of Medicolegal Death Investigators	
<u>Professional Services:</u>	Director of Laboratories Hunterdon Medical Center 2100 Wescott Dr. Flemington, NJ 08822	1/1999-present
	Medical Director LabCorp Clinical Trials Laboratory 750 Walnut Avenue Cranford, New Jersey 07016	7/2004-present
	Associate Director Miller-Keystone Memorial Blood Bank Bethlehem, PA 18017	1/1994-present
	Medical Examiner of Hunterdon County 2100 Wescott Drive Flemington, New Jersey 08822	1/2003-present
	Medical Director, Tandem Laboratories A Clinical Research Organization 115 Silva Street Ewing, New Jersey 08628	6/2008-present

Handwritten signature and date:
1/24/14

CURRICULUM VITAE
Steven Diamond, DO, MBA, FCAP

Publications:

1. Allen HB, Johnson BL, Diamond SM: Smoker's Wrinkles?: JAMA 225:1067
2. Diamond SM: Primary Adenocarcinoma of Renal Pelvis: Mil Med 139:959
3. Takeda M, Diamond SM: Cytological Diagnosis of Malignant Melanoma Metastatic to Endometrium: Acta Cytol 22:503

Employment History:

Warren Hospital
185 Roseberry Street
Phillipsburg, NJ 08865
7/1975-12/1998

Hunterdon Medical Center
2100 Wescott Drive
Flemington, NJ 08822
1/1999-present

SD
1/24/14

CV MEDICAL DIRECTOR – TANDEM LABS RTP

Home:
901 Arrowhead Lane
Mebane, NC 27302
Phone (919) 957-7446
E-mail: jbnofsinger@yahoo.com

Work:
Tandem Labs, Inc.
2202 Ellis Road
Durham, NC 27703
Phone (919) 281-0958

John Brian Nofsinger, Jr.

Experience

Laboratory Director


02/12 – Present Tandem Labs, Inc./LabCorp Durham, NC

Associate Laboratory Director

06/10 – 02/12 Tandem Labs, Inc./LabCorp Durham, NC

05/03 – 06/10 Enthalpy Analytical, Inc. Durham, NC

- Established GLP compliant bioanalytical group at Enthalpy Analytical, Inc.
- Developed LC/MS/MS methods for a wide array of molecules including peptides, chiral amino acids, drug interaction compounds, and an extended list of proprietary customer compounds
- Performed method and instrument troubleshooting activities
- Performed a gap analysis of Enthalpy's 21 CFR Parts 11 and 58 compliance. This analysis was then used to develop a road map toward regulatory compliance.
- Generated risk assessments, outlining the bioanalytical group's regulatory direction
- Developed and implemented GLP compliant system validations (DQ, IQ, and OQ) for Sciex API-4000 and API-5000 LC/MS/MS systems, Sciex's Analyst, Thermo's Watson laboratory information management system (LIMS), Zymark TurboVaps, Hewlett-Packard 8452A UV-Vis absorption spectrophotometers, a Perkin Elmer Multiprobe IIEX liquid handling robot, Hamilton Star liquid handling robots, a Veriteq computer controlled temperature monitoring and alarm system, a backup power system (Generac generator and Powerware UPS units), refrigerators, freezers, and various IBM and Dell personal computers
- Wrote and used standard operating procedures (SOPs) for cold storage unit and room temperature monitoring, software security configurations, building security, record keeping and notebooks, electronic signatures and audit trails, instrument use, instrument maintenance, a reference material program, data generation and use, reporting, study protocols, performance of phase inspections, master study schedule generation and maintenance, and data and record archival
- Formed the transition plan following the acquisition of Enthalpy Analytical's bioanalytical division by Tandem Labs/LabCorp
- Directly oversaw the analysis of hundreds of FDA regulated studies, including over 150 emergency IND studies for the monitoring of drug and metabolite concentrations aimed at tailoring patient medical care

 1313

Pharmaceutical Project Scientist/Project Scientist Coordinator

08/01 – 05/03

Triangle Laboratories, Inc.

Durham, NC

- Identified and solicited prospective clients while handling technical inquiries from current and prospective clients
- Developed new analytical chemistry assays using HPLC-UV and LC/MS/MS that met timely and strict regulatory standards
- Monitored GLP compliance in analytical projects
- Trained chemists and technicians in the application of new and modified analytical techniques

Education

Ph.D., Physical/Analytical Chemistry

08/97 – 08/01

Duke University

Durham, NC

- Graduate project aimed at correlating the photo-physics and photo-chemistry of melanins with molecular / microscopic structure
- Gained extensive experience in spectroscopies (steady-state absorption and emission, photoacoustic, time-correlated single photon counting, transient absorption), mass spectrometry (matrix-assisted laser desorption ionization), and microscopies (atomic force, confocal)

B.S., Chemistry with Physics Minor, *Cum Laude*

08/93 – 05/97

Wake Forest University

Winston-Salem, NC

- Undergraduate research project focused on characterizing real time secondary nucleation in chiral crystals
- Developed and constructed a continuous liquid sample introduction unit for a polarimeter

Affiliations

American Association of Pharmaceutical Scientists

American Society for Mass Spectrometry

American Physics Society

Phi Lambda Upsilon (National Chemistry Honor Society)

Sigma Pi Sigma (National Physics Honor Society)

Society of Physics Students

 1-3-13

Honors

Triangle Laboratories, Inc. Performance Award, 2002

William R. Krigbaum Chemistry Fellow, 2000

Duke University Conference Travel Fellow, 2000

Feature article in Science News, 156 (12): 190-191, September 18, 1999

Publications

- 1) Hopkins SC, Nofsinger JB, Allen MS, Koch P, Varney MA. In vivo saturation binding of GABA-A receptor ligands to estimate receptor occupancy using liquid chromatography/tandem mass spectrometry. *Biopharm. Drug Dispos.* 2009, 30, 9-20.
- 2) Liu Y, Kempf V, Nofsinger JB, Weinert EE, Rudnicki M, Wakamatsu K, Ito S, Simon JD. Comparison of the structural and physical properties of human hair eumelanin following enzymatic or acid/base extraction. *Pigment Cell Res.* 2003, 16, 355-365.
- 3) Nofsinger JB, Weinert EE, Simon JD. Towards the establishment of structure-function relationships for eumelanin. *Biospectroscopy.* **2002**, 67, 302-305.
- 4) Nofsinger JB, Liu Y, Simon JD. Aggregation of eumelanin mitigates photogeneration of reactive oxygen species. *Free Radical Bio. Med.* **2002**, 32, 720-730.
- 5) Nofsinger JB, Simon JD. Radiative relaxation of *Sepia* eumelanin is affected by aggregation. *Photochem. Photobiol.* **2001**, 74, 31-37.
- 6) Nofsinger JB, Ye T, Simon JD. Ultrafast nonradiative relaxation dynamics of eumelanin. *J. Phys. Chem. B.* **2001**, 105, 2864-2866.
- 7) Clancy CR, Nofsinger JB, Hanks RK, Simon JD. A hierarchical self-assembly of eumelanin. *J. Phys. Chem. B.* **2000**, 104, 7871-7873.
- 8) Nofsinger JB, Forest SE, Eibest LM, Gold KA, Simon JD. Probing the building blocks of eumelanins using scanning electron microscopy. *Pigment Cell Res.* **2000**, 13, 179-184.
- 9) Forest SE, Lam WC, Millar DP, Nofsinger JB, Simon JD. A model for the activated energy transfer within eumelanin aggregates. *J. Phys. Chem. B.* **2000**, 104, 811-814.
- 10) Nofsinger JB, Forest SE, Simon JD. Explanation for the disparity among absorption and action spectra of eumelanin. *J. Phys. Chem. B.* **1999**, 103, 11428-11432.

JS 1-3-13

Presentations

55th Annual Meeting for the American Society for Mass Spectrometry

Poster entitled *Determination of Epinephrine Concentrations in Human Plasma by High Performance Liquid Chromatography with Tandem Mass Spectrometry*. Nofsinger JB, Allen MS, Sykes DC, Hamzavi M. Indianapolis, IN (06/07).

29th Annual Meeting for the American Society for Photobiology

Talk entitled *Aggregation State-Dependent Eumelanin Photoreactivity*. Nofsinger JB, Clancy CR, Simon JD. Chicago, IL (07/01).

Symposium on Photoprotection

Talk entitled *Are There Structure-Function Relationships for Melanins?* Simon JD, Nofsinger JB, Clancy CR. Krakow, Poland (5/01).

2000 Duke University Biochemistry Retreat

Talk entitled *The Particle Size-Dependence of Eumelanin Photoreactivity*. Nofsinger JB, Clancy CR, Il'ichev YV, Simon JD.
Poster entitled *The Ultrafast, Primary Events in the Photoreaction of Eumelanin*. Ye T, Nofsinger JB, Simon JD. Wrightsville Beach, NC (09/00).

8th International Conference on Applications in Life Sciences

Talk entitled *Temporal and Spatial Studies of Melanins*. Simon JD, Nofsinger JB, Clancy CR, Il'ichev YV. Tokyo, Japan (08/00).

28th Annual Meeting for the American Society for Photobiology

Poster entitled *The Particle Size-Dependence of Eumelanin Photoreactivity*. Nofsinger JB, Il'ichev YV, Forest SE, Simon JD. San Francisco, CA (07/00).

XVII International Pigment Cell Conference

Posters entitled *The Size-Dependent Properties of Melanin and The Fluorescence Properties of Eumelanins*. Simon JD, Nofsinger JB, Forest SE. Nagoya, Japan (10/99).

XIX International Conference on Photochemistry

Poster entitled *The Size-Dependent Optical Properties of Melanins*. Nofsinger JB, Forest SE, Simon JD. Durham, NC (08/99).

27th Annual Meeting for the American Society for Photobiology

Talk entitled *Spectroscopic Studies of Size-Selected Melanin*. Nofsinger JB, Forest SE, Simon JD. Washington, DC (07/99).

60th Annual Meeting of the Society for Investigative Dermatology

Poster entitled *The Size-Dependent Optical Properties of Melanins*. Forest SE, Nofsinger JB, Simon JD. Chicago, IL (05/99).

113th North Carolina American Chemical Society Local Section Meeting

Talk entitled *Spectroscopic Studies of Size-Selected Melanin*. Simon JD, Nofsinger JB, Forest SE. Chapel Hill, NC (04/99).

43rd Annual Meeting of the Biophysical Society

Poster entitled *The Spectroscopy and Photoreactions of Natural and Synthetic Melanins*. Nofsinger JB, Forest SE, Studer-Martinez S, Drake KB, Simon JD. Baltimore, MD (02/99).

JB 1-3-13

50th Southeastern Regional Meeting of the American Chemical Society

Talk entitled *The Photophysics of Melanin*. Nofsinger JB, Forest SE, Studer-Martinez S, Drake KB, Simon JD. Research Triangle Park, NC (11/98).

8th Annual Meeting of the Panamerican Society for Pigment Cell Research

Poster entitled *Photochemistry of Melanin*. Forest SE, Nofsinger JB, Drake K, Simon JD. Snowmass, CO (08/98).

 1-3-13

CV MEDICAL DIRECTOR – TANDEM LABS NEW JERSEY

CURRICULUM VITAE

Steven Diamond, DO, MBA, FCAP

<u>Board Certification:</u>	American Board of Pathology Anatomic and Clinical Pathology (Recertified) Hematology	11/75 1/09 5/77
<u>Board Registry:</u>	America Board of Medicolegal Death Investigators	5/20/02
<u>Education:</u>	Chicago College of Osteopathic Medicine Midwestern University Chicago, IL	DO 6/70
	University of Wisconsin Madison, WI	Path Resident 1970 – 1971
	United States Navy Philadelphia, PA	Pathologist 1971 – 1973
	Thomas Jefferson University Philadelphia, PA	Path Resident 1973 – 1975
	Wharton School of Business University of Pennsylvania Philadelphia, PA	MBA 1992 – 1994
<u>State Licenses:</u>	New Jersey Pennsylvania Missouri New York North Carolina South Carolina Virginia Georgia Ohio Indiana	
<u>Medical Societies:</u>	College of American Pathologists Pennsylvania Medical Society National Association of Medical Examiners American Board of Medicolegal Death Investigators	
<u>Professional Services:</u>	Director of Laboratories Hunterdon Medical Center 2100 Wescott Dr. Flemington, NJ 08822	1/1999-present
	Medical Director LabCorp Clinical Trials Laboratory 750 Walnut Avenue Cranford, New Jersey 07016	7/2004-present
	Associate Director Miller-Keystone Memorial Blood Bank Bethlehem, PA 18017	1/1994-present
	Medical Examiner of Hunterdon County 2100 Wescott Drive Flemington, New Jersey 08822	1/2003-present
	Medical Director, Tandem Laboratories A Clinical Research Organization 115 Silva Street Ewing, New Jersey 08628	6/2008-present

1/24/14

CURRICULUM VITAE
Steven Diamond, DO, MBA, FCAP

Publications:

1. Allen HB, Johnson BL, Diamond SM: Smoker's Wrinkles?: JAMA 225:1067
2. Diamond SM: Primary Adenocarcinoma of Renal Pelvis: Mil Med 139:959
3. Takeda M, Diamond SM: Cytological Diagnosis of Malignant Melanoma Metastatic to Endometrium: Acta Cytol 22:503

Employment History:

Warren Hospital
185 Roseberry Street
Phillipsburg, NJ 08865

7/1975-12/1998

Hunterdon Medical Center
2100 Wescott Drive
Flemington, NJ 08822

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